



## Two-in-one suspension

DiaDuo-2 is an all-purpose, two-in-one suspension containing both diamonds and lubricant, for materialographic fine grinding and polishing.

You only have to apply DiaDuo-2 to the correct level of lubrication and you will never overdose the amount of diamonds. The relation between diamonds and lubricant is always the same regardless of dosage.

DiaDuo-2 is a stable suspension and it can be used problem-free in automatic dosing systems with peristaltic pumps.

**Diamond suspension and cooling lubricant combined into one product**

**Easy to clean after preparation**

**Does not contain solvents**

**Colour-coded**

**Contains monocrystalline diamonds for routine preparation**



## Easy to use

DiaDuo-2 is very easy to use for preparation and it is compounded to facilitate easy cleaning of the specimens after preparation. This makes DiaDuo-2 the best all-round diamond suspension available and the first choice, both for you who have non-complicated specimens and for you who have little or limited experience in materialographic preparation.

## Water-based suspension

DiaDuo-2 is water-based and does not contain solvents. The composition of DiaDuo-2 is harmless and does not require any labelling. Hence, DiaDuo-2 is an environmentally friendly choice.

## Colour-coded

DiaDuo-2 is colour-coded for easy identification of the diamond sizes. You will never mix up the different grain sizes.

## Contains monocrystalline diamonds

DiaDuo-2 contains monocrystalline diamonds which makes it the premier choice for routine preparation. It ensures reproducible results even when using manual dispensers, as the level between diamonds and lubricant is always correct.



### Specifications

Grain sizes	Contents	Code
1 µm	500 ml	40600036
	5 litre	40600041
3 µm	500 ml	40600037
	5 litre	40600042
6 µm	500 ml	40600038
	5 litre	40600043
9 µm	500 ml	40600039
	5 litre	40600044
15 µm	500 ml	40600040
	5 litre	40600045
<b>Manual dispensers</b>	<b>Volume</b>	<b>Code</b>
Dosing Gun	500 ml	40700031
Squeeze Bottle	250 ml	40700032



*Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.*



#### AUSTRALIA & NEW ZEALAND

**Struers Australia**  
18 Mayneview Street  
Milton QLD 4064  
Australia  
Phone: +61 7 3512 9600  
Fax: +61 7 3369 8200  
info.au@struers.dk

#### BELGIQUE (Wallonie)

**Struers S.A.S.**  
370, rue du Marché Rollay  
F-94507 Champigny  
sur Marne Cedex  
Téléphone +33 1 5509 1430  
Télécopie +33 1 5509 1449  
struers@struers.fr

#### CANADA

**Struers Ltd.**  
7275 West Credit Avenue  
Mississauga, Ontario L5N 5M9  
Phone +1 905-814-8855  
Fax +1 905-814-1440  
info@struers.com

#### CHINA

**Struers Ltd.**  
Rm. 702 Hi-Shanghai  
No. 970 Dalian Road  
Shanghai 200092, P.R. China  
Phone +86 (21) 5228 8811  
Fax +86 (21) 5228 8821  
struers.cn@struers.dk

#### CZECH REPUBLIC

**Struers GmbH**  
Organizační složka  
Havlíčková 361  
CZ-252 63 Roztoky u Prahy  
Phone +420 233 312 625  
Fax: +420 233 312 640  
czechrepublic@struers.de

#### DEUTSCHLAND

**Struers GmbH**  
Carl-Friedrich-Benz-Strasse 5  
D-47877 Willich  
Telefon +49 (0) 2154 486-0  
Fax +49 (0) 2154 486-222  
verkauf.struers@struers.de

#### FRANCE

**Struers S.A.S.**  
370, rue du Marché Rollay  
F-94507 Champigny  
sur Marne Cedex  
Téléphone +33 1 5509 1430  
Télécopie +33 1 5509 1449  
struers@struers.fr

#### IRELAND

**Struers Ltd.**  
Unit 11 Evolution@ AMP  
Whittle Way  
Catcliffe  
Rotherham S60 5BL  
Tel: 0845 604 6664  
Fax: 0845 604 6651  
info@struers.co.uk

#### ITALY

**Struers Italia**  
Via Monte Grappa 80/4  
20020 Arese (MI)  
Tel. +39-02/38236281  
Fax +39-02/38236274  
struers.it@struers.dk

#### JAPAN

**Marumoto Struers K.K.**  
Takara 3rd Building  
18-6, Higashi Ueno 1-chome  
Taito-ku, Tokyo 110-0015  
Phone +81 3 5688 2914  
Fax +81 3 5688 2927  
struers@struers.co.jp

#### NETHERLANDS/BELGIUM

**Struers GmbH Nederland**  
Electraweg 5  
NL-3144 CB Maassluis  
Tel. +31 (0) 10 599 72 09  
Fax +31 (0) 10 599 72 01  
netherlands@struers.de

#### POLAND

**Struers Sp. z o.o.**  
Oddział w Polsce  
ul. Jasnogórska 44  
31-358 Kraków  
Phone +48 12 661 20 60  
Fax +48 12 626 01 46  
poland@struers.de

#### ROMANIA

**Struers GmbH**  
Sucursala Sibiu  
Str.Scoala de Inot, nr. 18  
RO-550005 Sibiu  
Phone +40 269 244 558  
Fax +40 269 244 559  
romania@struers.de

#### SCHWEIZ

**Struers GmbH**  
Zweigniederlassung Schweiz  
Weissenbrunnenstrasse 41  
CH-8903 Birmsdorf  
Telefon +41 44 777 63 07  
Fax +41 44 777 63 09  
switzerland@struers.de

#### SINGAPORE

**Struers Singapore**  
627A Aljunied Road,  
#07-08 BizTech Centre  
Singapore 389842  
Phone +65 6299 2268  
Fax +65 6299 2661  
struers.sg@struers.dk

#### SUOMI

**Struers A/S, Sivuliike Suomessa**  
Kalevankatu 43  
00180 Helsinki  
Puhelin 0207 919 430  
Faksi 0207 919 431  
finland@struers.dk

#### SWEDEN

**Struers A/S**  
Ekbacksvägen 22, 3 tr  
SE-168 69 Bromma  
Telefon +46 (0)8 447 53 90  
Telefax +46 (0)8 447 53 99  
struers@struers.dk

#### UNITED KINGDOM

**Struers Ltd.**  
Unit 11 Evolution@ AMP  
Whittle Way  
Catcliffe  
Rotherham S60 5BL  
Tel: 0845 604 6664  
Fax: 0845 604 6651  
info@struers.co.uk

#### USA

**Struers Inc.**  
24766 Detroit Road  
Westlake, OH 44145-1598  
Phone +1 440 871 0071  
Fax +1 440 871 8188  
info@struers.com

#### ÖSTERREICH

**Struers GmbH**  
Zweigniederlassung Österreich  
Ginzkeyplatz 10  
A-5020 Salzburg  
Telefon +43 662 625 711  
Fax +43 662 625 711 78  
austria@struers.de

#### Struers A/S

Pederstrupvej 84  
DK-2750 Ballerup,  
Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk  
www.struers.com

