



The perfect match for faster results

- All-in-one suspension
- Can reduce preparation time by 30%
- Easy to use and clean




Faster preparation time and outstanding results

DiaPro is a line of unique diamond suspensions that make materialographic specimen preparation faster, more cost-efficient and easier to do.



Developed and optimised to deliver maximum performance


Extensive tests have shown that the composition of the surface, such as the polishing textile, has a remarkably high influence on the efficiency of the diamond suspensions used. Each DiaPro suspension has been developed and optimised to deliver maximum performance when used with Struers' specially designed grinding discs or polishing cloths. When used together, DiaPro suspensions and Struers surfaces give you a faster and more efficient way to realize reproducible and reliable results – every time.



Each DiaPro suspension is developed and optimised for different surface compositions.



Material removal chart



Reduce average preparation time by 30%

DiaPro suspensions are uniquely formulated to ensure a high removal rate. Tests show that DiaPro can reduce preparation time by 20-60% compared to traditional diamond products, depending on the used surface and the individual preparation method.


The higher removal rate means you can prepare specimens faster, and with less effort. You save time and consumables, and boost specimen preparation efficiency – without compromising quality.

7 Reasons to choose DiaPro


- High removal rate = Fast preparation time
- Consistent, reproducible results
- All-in-one
- No run-in time
- Easy cleaning
- Water-based, non-toxic
- Developed, produced and tested by Struers

Run-in time

When using traditional diamond products, the cloths usually need 5-15 min. run-in time, but with DiaPro suspensions, polishing cloths run-in time is practically eliminated.



Struers' cloths




Competing cloths



Make specimen preparation easier and more efficient

DiaPro is an easy to use all-in-one suspension meaning that the optimal mixing ratio of high-performance diamonds and lubricants has been selected carefully and combined into one product. You no longer need to worry about applying enough or overdosing the diamonds. Based on extensive research at Struers, the optimal ratio of diamonds and lubricant has already been calculated for you. Just apply DiaPro to achieve the correct level of lubrication.

High removal rate implies a lot of grinding and polishing dirt. DiaPro suspensions are composed to make it easier to clean the specimens, the surface and the machine.



* Specimen size \varnothing 30 mm

SPECIFICATIONS**Cat.no****DiaPro**

DiaPro are unique diamond suspensions for faster, more cost-efficient and much easier materialographic specimen preparation. Each single DiaPro suspension has been developed and optimized for use on specific grinding or polishing surfaces.

DiaPro Pan 15 µm

For pre-polishing on MD-Pan and DP-Pan

500 ml

5 l

40600375

40600376

DiaPro Allegro/Largo 9 µm

For fine grinding on MD-Allegro and MD-Largo

500 ml

5 l

40600369

40600370

DiaPro Plan 9 µm

For polishing on MD-Plan and DP-Plan

500 ml

5 l

40600385

40600386

DiaPro Sat 6 µm

For polishing on MD-Sat and DP-Sat

500 ml

5 l

40600377

40600378

DiaPro Largo 3 µm

For polishing on MD-Largo

500 ml

5 l

40600391

40600392

DiaPro Mol B 3 µm

For polishing of hard materials (>150 HV) on MD-Mol and DP-Mol

500 ml

5 l

40600379

40600380

DiaPro Mol R 3 µm

For polishing of soft materials (<150 HV) on MD-Mol and DP-Mol

500 ml

5 l

40600393

40600394

DiaPro Dur 3 µm

For polishing on MD-Dur and DP-Dur

500 ml

5 l

40600395

40600396

DiaPro Dac 3 µm

For polishing on MD-Dac and DP-Dac

500 ml

5 l

40600371

40600372

DiaPro Plus 3 µm

For polishing on MD-Plus and DP-Plus

500 ml

5 l

40600387

40600388

DiaPro Floc 3 µm

For polishing on MD-Floc and DP-Floc

500 ml

5 l

40600389

40600390

DiaPro Dur 1 µm

For polishing on MD-Dur and DP-Dur

500 ml

5 l

40600381

40600382

DiaPro Nap B 1 µm

For polishing of hard materials (>150 HV) on MD-Nap and DP-Nap

500 ml

5 l

40600373

40600374

DiaPro Nap R 1 µm

For polishing of soft materials (<150 HV) on MD-Nap and DP-Nap

500 ml

5 l

40600397

40600398

DiaPro Nap ¼ µm

For polishing on MD-Nap and DP-Nap

500 ml

5 l

40600383

40600384

Dosers

For manual dosing

Dosing Gun with 500 ml

Squeeze bottle 250 ml

40700031

40700032

For information on automatic dosing systems please require separate leaflet.

We care about the working environment

DiaPro suspensions are purely water-based and do not contain solvents. They are classified as non-toxic and thus the environmental-friendly choice of consumable. The products are extremely stable and do not settle, ensuring high quality and reproducibility.



The DiaPro suspension comes in 500 ml and 5 l.

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk
 www.struers.com



AUSTRALIA & NEW ZEALAND

Struers Australia
 27 Mayneview Street
 Milton QLD 4064
 Australia
 Phone +61 7 3512 9600
 Fax +61 7 3369 8200
 info.au@struers.dk

BELGIQUE (Wallonie)

Struers S.A.S.
 370, rue du Marché Rollay
 F- 94507 Champigny
 sur Marne Cedex
 Téléphone +33 1 5509 1430
 Télécopie +33 1 5509 1449
 struers@struers.fr

BELGIUM (Flanders)

Struers GmbH Nederland
 Elektraweg 5
 3144 CB Maassluis
 Telefoon +31 (10) 599 7209
 Fax +31 (10) 5997201
 netherlands@struers.de

CANADA

Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

CHINA

Struers Ltd.
 No. 1696 Zhang Heng Road
 Zhang Jiang Hi-Tech Park
 Shanghai 201203, P.R. China
 Phone +86 (21) 6035 3900
 Fax +86 (21) 6035 3999
 struers@struers.cn

CZECH REPUBLIC

Struers GmbH
 Organizační složka
 Havlíčkova 361
 CZ-252 63 Roztoky u Prahy
 Phone +420 233 312 625
 Fax +420 233 312 640
 czechrepublic@struers.de

DEUTSCHLAND

Struers GmbH
 Carl-Friedrich-Benz-Straße 5
 D-47877 Willich
 Telefon +49 (0) 2154 486-0
 Fax +49 (0) 2154 486-222
 verkauf@struers.de

FRANCE

Struers S.A.S.
 370, rue du Marché Rollay
 F-94507 Champigny
 sur Marne Cedex
 Téléphone +33 1 5509 1430
 Télécopie +33 1 5509 1449
 struers@struers.fr

IRELAND

Struers Ltd.
 Unit 11 Evolution@ AMP
 Whittle Way, Catcliffe
 Rotherham S60 5BL
 Tel. +44 0845 604 6664
 Fax +44 0845 604 6651
 info@struers.co.uk

ITALY

Struers Italia
 Via Monte Grappa 80/4
 20020 Arese (MI)
 Tel. +39-02/38236281
 Fax +39-02/38236274
 struers.it@struers.it

JAPAN

Marumoto Struers K.K.
 Takara 3rd Building
 18-6, Higashi Ueno 1-chome
 Taito-ku, Tokyo 110-0015
 Phone +81 3 5688 2914
 Fax +81 3 5688 2927
 struers@struers.co.jp

NETHERLANDS

Struers GmbH Nederland
 Elektraweg 5
 3144 CB Maassluis
 Telefoon +31 (10) 599 7209
 Fax +31 (10) 5997201
 netherlands@struers.de

ÖSTERREICH

Struers GmbH
 Zweigniederlassung Österreich
 Betriebsgebiet Puch Nord 8
 5412 Puch
 Telefon +43 6245 70567
 Fax +43 6245 70567-78
 austria@struers.de

POLAND

Struers Sp. z o.o.
 Oddział w Polsce
 ul. Jasnogórska 44
 31-358 Kraków
 Phone +48 12 661 20 60
 Fax +48 12 626 01 46
 poland@struers.de

ROMANIA

Struers GmbH
 Sucursala Sibiu
 Str.Scoala de Inot, nr. 18
 RO-550005 Sibiu
 Phone +40 269 244 558
 Fax +40 269 244 559
 romania@struers.de

SCHWEIZ

Struers GmbH
 Zweigniederlassung Schweiz
 Weissenbrunnstraße 41
 CH-8903 Birmensdorf
 Telefon +41 44 777 63 07
 Fax +41 44 777 63 09
 switzerland@struers.de

SINGAPORE

Struers Singapore
 627A Aljunied Road,
 #07-08 BizTech Centre
 Singapore 389842
 Phone +65 6299 2268
 Fax +65 6299 2661
 struers.sg@struers.dk

SUOMI

Struers Suomi
 Kalevankatu 43
 00180 Helsinki
 Puhelin +358 (0)207 919 430
 Faksi +358 (0)207 919 431
 finland@struers.fi

About Struers

Struers is the world leader in materialographic innovation. We develop, manufacture and distribute quality materialographic preparation solutions to help our customers realize maximum productivity and performance in laboratories and production facilities around the world. Struers is based in Copenhagen, Denmark with affiliates in 19 countries. We have employees globally, including experienced metallographers, qualified service personnel, and a worldwide network of dealers.

Learn more

Contact a Struers sales representative today or visit www.struers.com/diapro



e-Shop

Buy your consumables
 in Struers new e-Shop!
e-shop.struers.com